

SPECIFICATIONS		HI 8424
Range	рН	-2.00 to 16.00 pH
	mV	±699.9 mV; ±1999 mV
	Temperature	-20.0 to 120.0°C / -4.0 to 248.0°F
Resolution	рН	0.01
	mV	0.1 mV; 1 mV
	Temperature	0.1°C / 0.1°F
Accuracy (@20°C)	рН	±0.01
	mV	±0.2 mV; ±1 mV
	Temperature	±0.4°C/±0.8°F
Calibration	рН	automatic, one or two point with three memorized buffer values (pH 4.01, 7.01, 10.01)
	Offset	±1 pH
	Slope	75 to 110%
Input Impedance		10 ¹² Ohm
Temperature Compensation		automatic, -20 to 120°C (-4 to 248°F) or manual without temperature probe
pH Electrode		HI 1230B PEI body pH electrode with BNC connector and 1 m (3.3') cable (included)
Temperature Probe		HI 7662 stainless steel temperatures probe with 1 m cable (included)
Battery Type / Life		9V / approximately 150 hours of continuous use; auto-off after 20 minutes of non-use or disabled (user-selectable)
Environment		0 to 50°C (32 to 122°F); RH max 100%
Dimensions		164 x 76 x 45 mm (6.5 x 3.0 x 1.8")
Weight		180 g (6.3 oz.)

pH/ORP Meter

- Automatic calibration
- Automatic Temperature Compensation
- Automatic buffer recognition
- Three standard buffer values
- Calibration tags on LCD
- Low battery indicator
- Auto shut-off
- HOLD function
- °C and °F measurements

HI 8424 is one of the most popular pH meters on the market. This instrument is able to perform pH, ORP and temperature measurements with a high degree of accuracy and fast response.

Calibration is automatic at one or two points with three memorized buffer values (pH 4.01, pH 7.01 and pH 10.01). Once the instrument has been calibrated, the buffer values used during calibration are displayed with tags on the LCD. This feature keeps users informed of the current calibration and helps to avoid taking measurements that are out of range.

Exchange the pH probe for an ORP probe to obtain ORP readings in the mV range.

ORDERING INFORMATION

HI 8424 is supplied with HI 1230B pH electrode, HI 7662 temperature probe, HI 70004 pH 4.01 buffer solution sachet, HI 70007 pH 7.01 buffer solution sachet, HI 700661 electrode cleaning solution sachets (2), battery, protective case and instructions.

ELECTROD	ES	
HI 1230B	PEI body pH electrode with BNC connector and 1 m (3.3') cable	
HI 3230B	PEI body ORP electrode with platinum sensor, BNC connector and 1 m (3.3') cable	
HI 7662	Temperature probe	
SOLUTIONS		
HI 5004L HI 5007L HI 5010L HI 7061L HI 70300L	pH 4.01 buffer solution, 500 mL pH 7.01 buffer solution, 500 mL pH 10.01 buffer solution, 500 mL Electrode cleaning solution, 500 mL Electrode storage solution, 500 mL	

ACCESSORIES

HI 710015	Shockproof rubber boot, blue
HI 710016	Shockproof rubber boot, orange

NT (주)뉴텍계기 서울특별시 금천구 가산동 60-15 리더스타워 501호 02-868-8648(대) sales@yesnt.co.kr www.yesnt.co.kr

www.hannainst.com